



BlasTaq™ 2X PCR MasterMix

Cat. No. G895.5, G895

Store at -20°C.

Product Description

BlasTaq™ 2X PCR MasterMix is a ready-to-use MasterMix containing **abm's BlasTaq™** DNA Polymerase in a uniquely-formulated buffer with gel loading dye. This strategically-engineered, next generation Taq Polymerase provides **rapid extension rates and robust performance**. With specialized reaction conditions, this polymerase provides increased processivity, yields, and sensitivity, while shortening reaction times by up to 70%, compared to wild-type Taq DNA polymerase. BlasTaq™ has 5'-3' polymerase and 5'-3' exonuclease activities, lacks 3'-5' exonuclease activity, and produces 3'-dA-tailed amplicons. PCR products made with BlasTaq™ can be used with TA cloning vectors.

Product Component	Quantity	Cat. No.
BlasTaq™ 2X PCR MasterMix ¹	500 rxn (4 x 1.25 ml)	G895.5 ²
BlasTaq™ 2X PCR MasterMix ¹	1000 rxn (8 x 1.25 ml)	G895

¹ Buffer contains 1.5 mM Mg²⁺.

² G895.5 제품은 한국에서만 공급합니다.

* G895.5 제품은 한국에서만 공급합니다.

한국 대리점 : (주)이츠바이오

T : (02) 3462-8658 F : (02) 3462-8659

e : info@itsbio.co.kr www.itsbio.co.kr

Protocol

1. Mix individual components before use and assemble reaction on ice.

Component	Volume
2X BlasTaq™ PCR MasterMix	10.0 µl
Forward Primer (10 µM)	1 µl
Reverse Primer (10 µM)	1 µl
Template DNA	Variable (100 ng genomic DNA)
Nuclease-free H ₂ O	Up to 20 µl

2. Gently mix the reaction components and briefly centrifuge. Run thermocycling conditions for standard PCR:

Step	Temperature	Duration
Initial Denaturation ²	95°C	3 min
25 – 35 Cycles	95°C 60°C ³ 72°C	15 sec 15 sec 15 sec/kb
Final Extension	72°C	1 min

² For most applications, an initial 3 minute denaturation step at 95°C is sufficient. Increase to 5 minutes for high-GC or difficult templates.

³ BlasTaq™'s PCR buffer allows for primer annealing at 60°C for most primers and adjust only if needed.

3. After PCR, maintain the reaction at 4°C or store at -20°C until use.
4. Analyze the amplification products by agarose gel electrophoresis.
5. Visualize by ethidium bromide or SafeView™ (Cat No. **G108**) staining.

General Notes

- Specialized buffer for higher yields, sensitivity, and specificity compared to wild-type Taq polymerase.
- Decrease reaction times by 70% using specialized protocol.